



Premium*Life*

Overbed Table





INFECTION
CONTROL

Antimicrobial Thermoform top

PremiumLift antimicrobial table top slow the spread of microorganisms. Microorganisms include bacteria, viruses, protozoans, and fungi.

The thermoformed top is designed with infection control in mind. The surface is non-porous and gap free, which eliminates areas for water to collect and bacteria to grow.

PremiumLift Overbed Table unique design has an Auto-lift mechanism allowing the patient to effortlessly raise the table top from any position.

The table top will remain at the required position with no further adjustments and is lowered by using the handle with no bending required.



Split Tilt Top





Ground Floor

Cardiorespiratory

Endoscopy

Ward 2A

Ward 2B

Admissions

Accident and Emergency



Specifications

PremiumLift Overbed Table			
Table top	800mm x 400mm	Tilt Top Angle	0 – 70 Degrees
Split Tilt Top	600 - 200 x 400mm	Castors	50mm - 2 x Lockable
Low Height	780mm	Safe Working Load	20 kg
High Height	1050mm	Weight / Mass	15 kg

Functions & Features

STANDARD FEATURES

- Thermoform top
- Auto-Lift
- Lockable castors
- ISO Certified
- 3 year limited warranty

OPTIONAL FEATURES

- Raised edge & Split Tilt Top
- Color customization
- Low Height Frame



Premium *Lift*

Premium *Lift*

Address: 7/112 Colemans Road, Carrum Downs VIC 3201 Australia

T 1800 771 381

E sales@premiumlift.com.au

www.premiumlift.com.au

ARTG No: 310585

All parts and data continually undergo further development and may therefore deviate from the details given.

© Copyright PremiumLift
Coda Furniture Holding Pty Ltd

Our Full Terms & Conditions are available by request or can be found on our Website.
PremiumLift reserves the right to alter or amend this document without notice

Distributor Stamp:

