


PremiumLift

Bed Accessories




Self Help Pole

- Assist with moving and positioning in bed
- SWL Patient 80kg
- Powder coated steel frame and PVC handle
- Compatible with Ultra Low Bed and Low Low Bed **SKU: PLSHP**



King Single Extension Kit

- Works for both Low Low Bed and Ultra Low Bed **SKU: KSKIT**
- Quick and easy bolt on design
- The single bed frame can be easily transformed to a king single size without having to buy a new bed



Grab Rail

- Can be used as a support pole, assists the patient into a seated position
- 400mm x 400mm
- Screw securely onto the bed frame
- Strong powder coated steel frame
- Compatible with Ultra Low Bed and Low Low Bed (grey) **SKU: PLGR**
- Compatible with Deluxe Floor Bed (black) **SKU: PLDGR**



Length Extension

- 20cm Length Extension
- Screw securely onto the bed frame
- Strong powder coated steel
- Compatible with Ultra Low Bed (grey) **SKU: PLULBLE**
- Compatible with Deluxe Floor Bed (black) **SKU: PLDFBLE**



Folding Side Rails – Full Length

- Ensures user remains safely in the bed and prevents falls risk
- Quick release pin mechanism
- Screw securely onto the bed frame
- Strong powder coated steel frame
- Supplied in pair
- Compatible with Ultra Low Bed and Low Low Bed (grey) **SKU: PLFSR**
- Compatible with Deluxe Floor Bed (black) **SKU: PLDFBFSR**



Under Bed Light

The Underbed Light can be used in conjunction with Low Low Bed **SKU: LLBUBL** and Deluxe Floor Bed **SKU: DFBUBL**. Once movement has been detected, a soft glow will illuminate around the bed side.

Backup Battery

The PremiumLift Backup Battery is compatible with our PremiumLift Bed range. In the event of a power failure, the backup battery will provide power. **SKU: PLBB**



Remote Control Pocket Holder

Compatible with our PremiumLift Deluxe Floor Bed **SKU: PLDFBH**



Upholstered Head & Foot Boards

- Ultra Low Bed Vanilla (beige) head & foot board set **SKU: PLULBBFB**
- Ultra Low Bed Classic French Grey (grey) head & foot board set **SKU: PLULBGFB**
- Deluxe Floor Bed Classic French Grey (grey) head board only **SKU: PLDFBHB**
- Deluxe Floor Bed Classic French Grey (grey) foot board only **SKU: PLDFBFB**

