


# Premium *Lift*

## Bedside Locker



ESSENTIAL SOLUTION

**User Friendly**



## **Economical and Reliable**

PremiumLift Bedside Locker design is not only a matter of appearance but also it affects the capacity of offering user friendly and functional product. This bedside cabinet features one top lockable drawer providing secure storage for medication and/or for valuable personal belongings. This cabinet is strong, durable and fitted with four castors.

PremiumLift Bedside Locker is designed with infection control in mind, the surface is easy to clean and moisture resistant allowing use of most harsh cleaning chemicals. The overall design is modern and will fit in well with any room decor.

**PremiumLift**





## Available Custom Design



**TGA** Approved Medical Device

## Technical Specification

✓ Bedside Locker SKU	PLBL303H
✓ Height	820mm
✓ Width	450mm
✓ Depth	500mm
✓ 3 Drawers	Yes
✓ Lock	Keyed Alike
✓ Board Thickness	16mm
✓ Colour	Timber Look Beech
✓ Drawer SLW	5kg (On maximum load only open one drawer at a time)
✓ Unit weight	31kg
✓ Castors	Ø 50mm
✓ Handles	Round Aluminium
✓ ISO Certified	Yes
✓ TGA Approved	Yes
✓ Warranty	3 year limited
✓ Packaging Size	90 x 50 x 55cm
✓ Packaging Weight	32kg



# Premium *Lift*



# Premium *Lift*

**Address:** 1 Force Way Pakenham VIC 3810 Australia

T 1800 771 381

E [sales@premiumlift.com.au](mailto:sales@premiumlift.com.au)

[www.premiumlift.com.au](http://www.premiumlift.com.au)

All parts and data continually undergo further development and may therefore deviate from the details given.

© **Copyright PremiumLift**

Our Full Terms & Conditions are available by request or can be found on our website. PremiumLift reserves the right to alter or amend this document without notice

Distributor Stamp:

