


Premium *Lift* Catalogue



Floorline Beds Specialists

Premium *Lift*

Deluxe Floor Bed



18cm from the floor for complete user and patient safety



Australia Lowest Home Care Bed

We have spoken with hundreds of Occupational Therapists, carers and aged care facility managers who express their need for a modern, cosy fully functional bed with floor-line capabilities and ergonomic height.

We were asked to design a bed with a hidden castors system in order to offer the ultimate floorline care bed that looks like a home bed.

We quickly discovered that carers and users alike are tired of buying substandard beds with limited functions, clinical look and unsafe design. They told us that poorly designed beds can influence mood, behaviour and most importantly the desire to interact with friends or family.

In 2022, together with the world leading supplier of electric bed components, we were able to design a revolutionary bed frame that can be lowered to just 18cm from the floor level and can elevate to an incredible height of 73cm and meet OH&S height standard requirements.

Additionally, we successfully developed a unique Hidden Castors System that will revolutionise the home care market and truly enhances patient's recovery whilst promoting emotional comfort.

Deluxe Floor Bed modern 'column-less' design and unique features offers 12 pre-set positions and over 60 position combinations, including the popular Zero Gravity and Trendelenburg functions. This bed will promote the carer's wellbeing, prevent injury and offer a memorable patient experiences and life-changing care.



Lowest Home Care Bed



Wireless Remote Control



Elegant Design









Hidden Castors System



12 Preset Positions



Cosy & Fully Functional

Standard Features & Advantages	Uses & Benefits Lowest Home Care Bed	Options Accessories & Extras
<ul style="list-style-type: none"> ✓ Min Height 18cm ✓ Hidden Castors System ✓ Non-Clinical Design ✓ 12 pre-set functions ✓ TGA Approved ✓ ISO Certified ✓ 10 year limited warranty ✓ Zero Gravity Function ✓ Seating Position ✓ 9V Battery Backup ✓ Elegant Head & Foot Boards ✓ SWL 200kg 	<div style="text-align: center;">  Fall Prevention </div> <div style="text-align: center;">  After Stroke Care </div> <div style="text-align: center;">  Alzheimer Care </div> <div style="text-align: center;">  Disability Care </div> <div style="text-align: center;">  Carers Support </div> <div style="text-align: center;">  Injury Recovery </div>	<ul style="list-style-type: none"> ✓ Under Bed Light ✓ 20cm Length Extension ✓ Remote Control Holder ✓ Foldable Side Rails ✓ Grab Rails ✓ Head & Foot Boards ✓ Pillow ✓ Fall Mat ✓ Mattress ✓ External Battery Pack





Australia Lowest Home Care Bed

18cm from the floor for complete user and patient safety



Hidden Castors System

Developed to offer the ultimate floor-line home care bed that looks like a home bed



TGA Approved Medical Device

Technical Specification

✓ King Single Size SKU	PLDFBK
✓ Overall Length	2150mm
✓ Height Range	180mm – 730mm
✓ King Single Dimensions	2050mm x 1070mm
✓ Length with Extension	2250mm
✓ Electric backrest	0° to 70°
✓ Electric legrest	0° to 38°
✓ Trendelenburg	20°
✓ Zero Gravity	Yes
✓ Seating Position	Yes
✓ Wireless Remote	Yes
✓ Under Bed Light	Yes
✓ Safe working load (SWL)	200kg
✓ Bed weight	80kg
✓ Hidden Castors	Ø 50 mm
✓ Electrics Brand	OKIN
✓ Battery Backup	Yes (9V)
✓ ISO Certified	Yes
✓ TGA Approved	Yes
✓ Warranty	10 year limited

Premium *Lift*

Ultra Low Bed



6.7cm from the floor for complete user and patient safety



Australia Lowest Aged Care Bed

PremiumLift Ultra Low Bed is the ultimate floor-line bed with a minimum height of just 6.7cm from the floor ensuring complete user and patient safety, suitable for residential aged care facilities and assisted living solutions.

Ultra Low Bed is designed in Germany to the higher quality specification and equipped with a state of the art Anti-Entrapment System. This safety mechanism will automatically reduce the lowering speed when the bed reaches 120mm from floor level, to allow the operation to be stopped and the bed raised before any limb or object becomes trapped, ensuring complete user and patient safety.

PremiumLift Ultra Low Bed have plenty of user-friendly accessories and boards that can be easily fitted to the bed frame offering the ultimate customer experience, support carers with everyday tasks and bring a new level of confidence and dignity to the user.









Timber Look Beech



Classic French Grey



Vanilla Beige

Standard Features & Advantages	Uses & Benefits Australia Lowest Bed	Options Accessories & Extras
<ul style="list-style-type: none"> ✓ Min Height 6.7cm ✓ Anti-Entrapment System ✓ Five-function electric bed ✓ TGA Approved ✓ ISO Certified ✓ 10 year limited warranty ✓ 4 section bed for easy assembly and dismantle ✓ 9V Battery Backup ✓ Timber Look Beech Head & Foot Boards ✓ SWL 185kg 	<div style="text-align: center;">  Fall Prevention </div> <div style="text-align: center;">  After Stroke Care </div> <div style="text-align: center;">  Alzheimer Care </div> <div style="text-align: center;">  Disability Care </div> <div style="text-align: center;">  Carers Support </div> <div style="text-align: center;">  Injury Recovery </div>	<ul style="list-style-type: none"> ✓ King Single Width Extension Kit ✓ 20cm Length Extension ✓ Fabric Head and Foot Boards ✓ Foldable Side Rails ✓ Grab Rail ✓ Self Help Pole ✓ Fall Mat ✓ Mattress ✓ External Battery Pack





Australia Lowest Bed

6.7cm from the floor for complete user and patient safety



King Single Width Extension

Single to King Single within in seconds



TGA Approved Medical Device

Technical Specification

✓ Single Size SKU	PLULB
✓ King Single Size SKU	PLULBK
✓ Overall Length	2150mm
✓ Height Range	67mm – 650mm
✓ Single Dimensions	2000mm x 900mm
✓ King Single Dimensions	2000mm x 1070mm
✓ Length with Extension	2200mm
✓ Electric backrest	0° to 70°
✓ Electric legrest	0° to 38°
✓ Foot rest	0° to 30° (4 stages)
✓ Trendelenburg	20°
✓ Reverse Trendelenburg	20°
✓ Safe working load (SWL)	185kg
✓ Bed weight	80kg
✓ Locking Castors	Ø 75 mm
✓ Electrics Brand	Limoss
✓ Battery Backup	Yes (9V)
✓ ISO Certified	Yes
✓ TGA Approved	Yes
✓ Warranty	10 year limited

Premium *Lift*

Low Low Bed



20cm from the floor for complete user and patient safety




Reliable Low Profile Care Bed

PremiumLift Low Low Bed is designed to encourage the healing process by reducing discomfort when moving from the bed and preventing bed-fall injuries. It offers a variety of standard functions including 20cm low height, Trendelenburg, cardiac chair position and built-in 20cm length extension for excellent care to be provided.

We understand that recovery is extremely important among the elderly. This bed will allow the person you love to rest with comfort and security, preventing any potential injury from falling out of bed. PremiumLift Low Low Bed folding side rails will ensure that the user remains safely in bed.

PremiumLift Ultra Low Bed have plenty of user-friendly accessories that can be easily fitted to the bed frame offering the ultimate customer experience, support carers with everyday tasks and bring a new level of confidence and dignity to the user.









Reliable & Safe Design



User Friendly Accessories



8 Preset Positions

Standard Features & Advantages	Uses & Benefits Australia Lowest Bed	Options Accessories & Extras
<ul style="list-style-type: none"> ✓ Min Height 20cm ✓ Built-In Length Extension ✓ Five-function electric bed ✓ TGA Approved ✓ ISO Certified ✓ 10 year limited warranty ✓ 4 section bed for easy assembly and dismantle ✓ 9V Battery Backup ✓ Timber Look Beech Head & Foot Boards ✓ SWL 185kg 	<div style="text-align: center;">  Fall Prevention </div> <div style="text-align: center;">  After Stroke Care </div> <div style="text-align: center;">  Alzheimer Care </div> <div style="text-align: center;">  Disability Care </div> <div style="text-align: center;">  Carers Support </div> <div style="text-align: center;">  Injury Recovery </div>	<ul style="list-style-type: none"> ✓ King Single Width Extension Kit ✓ Foldable Side Rails ✓ Grab Rail ✓ Self Help Pole ✓ Fall Mat ✓ Mattress ✓ External Battery Pack






Reliable and Safe Design

20cm from the floor for complete user and patient safety



King Single Width Extension

Single to King Single within in seconds



TGA Approved Medical Device

Technical Specification

✓ Single Size SKU	PLLLB
✓ King Single Size SKU	PLLLBK
✓ Overall Length	2150mm
✓ Height Range	200mm – 650mm
✓ Single Dimensions	1950mm x 900mm
✓ King Single Dimensions	1950mm x 1070mm
✓ Length with Extension	2150mm
✓ Electric backrest	0° to 70°
✓ Electric legrest	0° to 38°
✓ Foot rest	0° to 30° (4 stages)
✓ Trendelenburg	20°
✓ Reverse Trendelenburg	20°
✓ Safe working load (SWL)	185kg
✓ Bed weight	80kg
✓ Locking Castors	Ø 75 mm
✓ Electrics Brand	Limoss
✓ Battery Backup	Yes (9V)
✓ ISO Certified	Yes
✓ TGA Approved	Yes
✓ Warranty	10 year limited

PremiumLift

Bed Accessories




Self Help Pole

- Assist with moving and positioning in bed
- SWL Patient 80kg
- Powder coated steel frame and PVC handle
- Compatible with Ultra Low Bed and Low Low Bed **SKU: PLSHP**



King Single Extension Kit

- Works for both Low Low Bed and Ultra Low Bed **SKU: KSKIT**
- Quick and easy bolt on design
- The single bed frame can be easily transformed to a king single size without having to buy a new bed



Grab Rail

- Can be used as a support pole, assists the patient into a seated position
- 400mm x 400mm
- Screw securely onto the bed frame
- Strong powder coated steel frame
- Compatible with Ultra Low Bed and Low Low Bed (grey) **SKU: PLGR**
- Compatible with Deluxe Floor Bed (black) **SKU: PLDGR**



Length Extension

- 20cm Length Extension
- Screw securely onto the bed frame
- Strong powder coated steel
- Compatible with Ultra Low Bed (grey) **SKU: PLUBLE**
- Compatible with Deluxe Floor Bed (black) **SKU: PLDFBLE**



Folding Side Rails – Full Length

- Ensures user remains safely in the bed and prevents falls risk
- Quick release pin mechanism
- Screw securely onto the bed frame
- Strong powder coated steel frame
- Supplied in pair
- Compatible with Ultra Low Bed and Low Low Bed (grey) **SKU: PLFSR**
- Compatible with Deluxe Floor Bed (black) **SKU: PLDFLR**

Under Bed Light

The Underbed Light can be used in conjunction with Low Low Bed **SKU: LLUBL** and Deluxe Floor Bed **SKU: DFBUBL**. Once movement has been detected, a soft glow will illuminate around the bed side.



Backup Battery

The PremiumLift Backup Battery is compatible with our PremiumLift Bed range. In the event of a power failure, the backup battery will provide power. **SKU: PLBB**

Remote Control Pocket Holder

Compatible with our PremiumLift Deluxe Floor Bed **SKU: PLDFBRCPH**

Upholstered Head & Foot Boards

- Ultra Low Bed Vanilla (beige) head & foot board set **SKU: PLULBBFB**
- Ultra Low Bed Classic French Grey (grey) head & foot board set **SKU: PLULBGFB**
- Deluxe Floor Bed Classic French Grey (grey) head board only **SKU: PLDFBHB**
- Deluxe Floor Bed Classic French Grey (grey) foot board only **SKU: PLDFBFB**




PremiumLift

Foam Mattresses

Multi-Flex

PremiumLift Multi-Flex mattress is designed with a multi-aeriation channel system that adapt to the user body weight, stimulate circulation through massage points and prevents bedsores.



Applications	Stimulate Circulation
Layers	Bottom Firm, Top Soft
Dimensions	90x195cm & 106x200cm
Thickness	15cm
Single Size SKU	PLMFS
K. Single Size SKU	PLMFKS
Max User Weight	150 kg
Cover Material	Four-way stretch removable, waterproof, fire retardant cover and foam



Premium-Flex

PremiumLift Premium-Flex mattress is designed with a three-layer premium foam construction that adapt to the user body weight.



Applications	Low to Medium Risk
Layers	3 Layers
Dimensions	90x195cm & 106x200cm
Thickness	15cm
Single Size SKU	PLPFS
K. Single Size SKU	PLPFKS
Max User Weight	200 kg
Cover Material	Four-way stretch removable, waterproof, fire retardant cover and foam



Bari-Flex

PremiumLift Bari-Flex mattress is designed with a three-layer premium foam construction that adapt to the user body weight.

Applications	Low to Medium Risk
Layers	3 Layers
Dimensions	90x195cm & 106x200cm
Thickness	15cm
Single Size SKU	PLBARFS
K. Single Size SKU	PLBARFKS
Max User Weight	350 kg
Cover Material	Four-way stretch removable, waterproof, fire retardant cover and foam



Fall Mat

PremiumLift Fall Mat helps reduce impact forces by up to 90% and provide full coverage of the bed exit area. Tri-fold design for easy storage

Applications	Injury Prevention
Layers	1 Layer
Dimensions	90x178cm
Thickness	3cm
SKU	PLFM
Max User Weight	200 kg
Cover Material	Tri-Fold, four-way stretch waterproof, fire retardant cover and foam




Home Care Mattress



Deluxe-Flex

PremiumLift Deluxe-Flex Home Care mattress is designed with a two layers premium foam construction and top gel memory foam layer that adapt to the user body weight. Suitable for Deluxe Floor Bed

Applications	Home Care
Layers	3 Layers
Dimensions	106x200cm
Thickness	20cm
K. Single Size SKU	PLDFKS
Max User Weight	200 kg
Cover Material	Hypoallergenic, fire retardant cover and foam



AirCare

Alternating Air Mattresses

AirCare Max-8 Alternating Air Mattress

The AirCare Max-8 Alternating Air Mattress is a high-quality replacement mattress. Designed for a frail user who could be vulnerable to developing pressure ulcers, the system provides a very high degree of pressure relief whilst protecting the user. It is recommended for use by a patient who is at very high-risk of developing pressure ulcers, and support the treatment of up to stage 4 pressure ulcers.

Features & Advantages

- ✓ Advanced Ultra-Quiet Pump System
- ✓ User-Friendly Smart Operation Panel
- ✓ 18 individual cells, 15 alternating and 3 static pillow section
- ✓ Individually replaceable air cells
- ✓ Static/ Max Inflate for Nurse use
- ✓ Adjustable user comfort level
- ✓ Safety digital pressure level indicator
- ✓ Low pressure and power failure alarm
- ✓ Advanced pressure detection technology
- ✓ Quick connect and release plug
- ✓ Taiwan made Nylon TPU cover 2-way stretch
- ✓ 50mm integral foam base for extra support
- ✓ Rapid release CPR valve for quick deflating
- ✓ Pump Wireless Remote Control

Maximum Patient Weights:	210kg
Single Dimensions:	900 x 1980 x 200mm (SKU ACMAX8S)
K. Single Dimensions:	1060 x 1980 x 200mm (SKU ACMAX8KS)
Mattress Weight:	7kg
Warranty:	2 Years
Pump:	AirCare Pro (SKU ACPP)
Pump Dimensions:	345 x 115 x 185mm




AirCare Max-8 Alternating Air Mattress - Replacement System



Rapid release CPR valve




18 individual cells



Quick connect and release plug



50mm integral foam base



AirCare Pro-Pum - Control Panel

Premium *Lift*

Overbed Table



ERGONOMIC AUTO-LIFT DESIGN

NO BENDING REQUIRED



Ergonomic Design & Antimicrobial Top

PremiumLift Overbed Table antimicrobial table top slow the spread of microorganisms. Microorganisms include bacteria, viruses, protozoans, and fungi.

The thermoformed top is designed with infection control in mind. The surface is non-porous and gap free, which eliminates areas for water to collect and bacteria to grow.

PremiumLift Overbed Table unique design has an Auto-lift mechanism allowing the patient to effortlessly raise the table top from any position.

The table top will remain at the required position with no further adjustments and is lowered by using the handle with no bending required.








Timber Look (Beech)



Polar White (White)



Pearl Grey (Light Grey)

Standard Features & Advantages	Uses & Benefits Australia Lowest Bed	Options Accessories & Extras
<ul style="list-style-type: none"> ✓ Ergonomic Design ✓ Auto-Lift Mechanism ✓ No Bending Required ✓ TGA Approved ✓ ISO Certified ✓ 3 year limited warranty ✓ Antimicrobial Top ✓ Split Tilt Top ✓ Timber Look or White Top ✓ SWL 20kg 	<div style="text-align: center;">  <p>After Stroke Care</p> </div> <div style="text-align: center;">  <p>Alzheimer Care</p> </div> <div style="text-align: center;">  <p>Disability Care</p> </div> <div style="text-align: center;">  <p>Carers Support</p> </div> <div style="text-align: center;">  <p>Injury Recovery</p> </div>	<ul style="list-style-type: none"> ✓ Raised Edge Table Top ✓ Table Top Colour ✓ Table Top Size ✓ Table Top Material ✓ Castors





Auto-Lift Mechanism

The unique design allows the user to adjust the height with one hand no bending is required



Infection Control

Antimicrobial table top slow the spread of microorganisms (Gap Free)



Technical Specification

✓ Non Tilt Top SKU	PLFOBT
✓ Split Tilt Top SKU	PLOBT
✓ Overall Length	800mm
✓ Height Range	750mm – 1000mm
✓ Table Top Dimensions	400mm x 800mm
✓ Tilt Top Fix Section	400mm x 200mm
✓ Tilt Top Tilt Section	400mm x 600mm
✓ Tilt Angle	0° to 70°
✓ Frame Colour	Grey
✓ Safe working load (SWL)	20kg
✓ Table weight	15kg
✓ Castors (2x lockable)	Ø 50mm nylon
✓ Lifting Mechanism	Gas Strut
✓ ISO Certified	Yes
✓ TGA Approved	Yes
✓ Warranty	3 year limited
✓ Packaging Size	83 x 52 x 15cm
✓ Packaging Weight	15.5kg

TGA Approved Medical Device


Premium *Lift*

Bedside Locker



ESSENTIAL SOLUTION

User Friendly



Economical and Reliable

PremiumLift Bedside Locker design is not only a matter of appearance but also it affects the capacity of offering user friendly and functional product. This bedside cabinet features one top lockable drawer providing secure storage for medication and/or for valuable personal belongings. This cabinet is strong, durable and fitted with four castors.

PremiumLift Bedside Locker is designed with infection control in mind, the surface is easy to clean and moisture resistant allowing use of most harsh cleaning chemicals. The overall design is modern and will fit in well with any room decor.

PremiumLift



Available Custom Design



TGA Approved Medical Device

Technical Specification

✓ Bedside Locker SKU	PLBL303H
✓ Height	820mm
✓ Width	450mm
✓ Depth	500mm
✓ 3 Drawers	Yes
✓ Lock	Keyed Alike
✓ Board Thickness	16mm
✓ Colour	Timber Look Beech
✓ Drawer SLW	5kg (On maximum load only open one drawer at a time)
✓ Unit weight	31kg
✓ Castors	Ø 50mm
✓ Handles	Round Aluminium
✓ ISO Certified	Yes
✓ TGA Approved	Yes
✓ Warranty	3 year limited
✓ Packaging Size	90 x 50 x 55cm
✓ Packaging Weight	32kg



Premium *Lift*



Premium *Lift*

Address: 1 Force Way Pakenham VIC 3810 Australia

T 1800 771 381

E sales@premiumlift.com.au

www.premiumlift.com.au

All parts and data continually undergo further development and may therefore deviate from the details given.

© **Copyright PremiumLift**

Our Full Terms & Conditions are available by request or can be found on our website. PremiumLift reserves the right to alter or amend this document without notice

Distributor Stamp:

